



Multiple-Operator hood

HI4P6S

Handling of chemicals

Convenient for heavy solvents or light smokes



Technical characteristics	
Débit Variable	50 to 450 m ³ /h
Capacity	72 W
L x H x W mm	430 x 350 x 200
Weight	5.7 kg
Power supply	230 V / 50 Hz
Maintenance	Periodic change of the filter
Exit evacuation	Yes, Ø125mm

Strengths

- Portable
- 6 extraction inputs
- Different accessories
- Delivered with 6 x 3m of ringed PVC pipe

Utilisation

Recycling : Active carbon filtration or dust filter

External output : external output filter + exhaust duct Ø 125 mm

Options



Suction plate AP1

- ∞ Pas d'installation à faire, se pose directement sur le poste de travail
- ∞ connects simply with a Ø50 mm pipe
- ∞ Dimensions 410 x 350 x 200 mm
- ∞ The adjustable visor helps channel flows



Flexible extraction arms

- ∞ For mobile or temporary installation
- ∞ Manually adjustable
- ∞ Different arrangements (on table, wall, or ceiling)
- ∞ Possibility to choose the materials : Standard, ESD, ATEX, polypropylen...)