



Overhead fume hood

HI9P

Solvents treatment
Handling of chemicals
Captation of light fumes

Technical characteristics

Débit Variable	50 to 450 m ³ /h
Capacity	80 W
L x H x W mm	720 x 640 x 300
Weight	8 kg
Power supply	230 V / 50 Hz
Maintenance	Periodic change of the filter
Exit evacuation	Yes, Ø125 mm

Strengths

Workstation under the fume hood
 Integrated lighting
 Excellent visibility

Utilisation

Recycling : Active carbon filtration or dust filter

External output : external output filter + exhaust duct Ø 125 mm

