

## Suitcase hood – HI3P

 Mobile extractor for handling chemicals, solvents, soldering fumes, small particles, glues etc.




- ✓ Portable fume hood
- ✓ Important suction surface
- ✓ Does not encumber the working plan
- ✓ Integrated lighting
- ✓ Can be equipped with visors and transformed into small cabinets

### Technical data

Variable flow	50 to 450 m <sup>3</sup> /h
Power	80 W
L x H x D	350 x 430 x 200 mm
Weight	6 kg
Mains supply	230 V / 50 Hz
Maintenance	Periodic change of the filter
Outer rejection	Yes, Ø125 mm

HI3P can be used both ways :

-  As a recirculating (ductless) fume hood with an active carbon filter or a particle filter
  - Or as an external rejection system using a dust filter and an exhaust duct Ø125 mm
- Possibility of using protection prefilters to lengthen the life of activated carbon filters

### Accessories and consumables



Small visor VHI3PP



Funnel visor VHI3GP



Active carbon filter FCP1

How to choose a filter ?

Find all our filters on [www.ftm-technologies.com](http://www.ftm-technologies.com)