

Suitcase hood – HI3P

 Mobile extractor for handling chemicals, solvents, soldering fumes, small particles, glues etc.




- ✓ Portable fume hood
- ✓ Important suction surface
- ✓ Does not encumber the working plan
- ✓ Integrated lighting
- ✓ Can be equiped with visors and transformed into small cabinets

Technical data

Variable flow	50 to 450 m3/h
Power	80 W
L x H x D	350 x 430 x 200 mm
Weight	6 kg
Mains supply	230 V / 50 Hz
Maintenance	Periodic change of the filter
Outer rejection	Yes, Ø125 mm

HI3P can be used both ways :

-  As a recirculating (ductless) fume hood with an active carbon filter or a particle filter
Or as an external rejection system using a dust filter and an exhaust duct Ø125 mm
Possibility of using protection prefilters to lengthen the life of activated carbon filters

Accessories and consumables



Small visor VHI3PP



Funnel visor VHI3GP



Active carbon filter FCP1

How to choose a filter ?

Find all our filters on www.ftm-technologies.com